

RegenaVate[®] DBM

“Flowable” DBM provides unique handling



REGENERATIVE SOLUTIONS

1 Exceptional Handling

- RegenaVate DBM Fill is “flowable” at 45°C and is a resilient solid at body temperature
- Will not wash away during irrigation
- RegenaVate DBM Fill allows the clinician to control product consistency: gel, paste or putty

3 Convenient & Easy To Use

- Packaged in syringes for convenient dispensing
- Easy rehydration

2 Induces Bone Formation & Facilitates Bone Growth*

- RegenaVate DBM Fill contains human demineralized bone matrix in an inert porcine gelatin carrier
- Contains 24 percent DBM by weight
- DBM is tested for osteoinductivity in a scientifically-proven in vivo rat assay

* These implants were evaluated in a human clinical study and were shown to induce bone formation. Each lot is tested using the athymic rat assay to verify osteoinductive potential.



ZIMMER BIOMET
Your progress. Our promise.®

Simple Handling and Regeneration



INDICATIONS

RegenaVate DBM Fill is intended to be packed into bony voids or gaps to fill and/or augment dental intraosseous, oral and cranio-/maxillofacial defects including:

- Extraction Sockets (implant preparation/placement)
- Alveolar Ridge Augmentation
- Sinus Floor Elevation

Ordering Information:

Item Number	Description
013005Z	RegenaVate DBM Fill RT, 0.5 cc
013010Z	RegenaVate DBM Fill RT, 1 cc

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. RegenaVate DBM is manufactured by RTI Surgical, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0900 REV A 10/19 ©2019 Zimmer Biomet. All rights reserved.

