

Quality & Confidence – ALL IN ONE



The **Screw Fixation System** provides a compact solution for the temporary fixation and stabilization of bone transplants, suitable resorbable and non-resorbable bone replacement materials, and membranes for the alveolar ridge. Two color coded systems in 1.5mmD MICRO (blue) and 2.0mmD MINI (red) make easy and rapid identification of the parts possible and simplifies parts matching. This modular storage system permits individual configuration and the open design aims to ensure convenient access during cleaning and sterilization. Fixation screws and mesh are manufactured from pure titanium or titanium alloy. They are biocompatible, corrosion-proof and nontoxic in the biological environment. They allow imaging virtually free of artifacts.



SECURE & SIMPLE

The **Screw Fixation System** comes with power grip for secure and stable transfer to the surgical site. The screws are easily picked up and have a reliable connection to the driver.

COMPACT

The screw cartridges of the **Screw Fixation System** are more than simply storage containers. The cartridges make organization and sterilization of screws easy. An imprinted four data identification system provides traceability.

SCREW FIXATION KIT – START-UP KIT ORDERING INFORMATION




Contents	Tray	Screw Driver Handle	Screw Driver Insert, Short	Screw Driver Insert, Long	Pilot Drill, Micro, 14mmL	Pilot Block Drill, Micro
Assembled Start-Up Kit*	69.01.10Z					
Components	69.01.11Z	75.23.52Z	75.23.23Z	75.23.19Z	69.01.09Z	69.01.16Z

*Screws are not included in the start-up kit



1.5mmD MICRO System	System/Color Indication	Color Marking	Measurements	Kit & Stand Alone	Stand Alone Only
Tray	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	69.01.11Z	
Screw Driver Handle	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	75.23.52Z	
Screw Driver Insert, Short	1.5 Micro; Blue	blue	17mmL	75.23.23Z	
Screw Driver Insert, Long	1.5 Micro; Blue	blue	55mmL	75.23.19Z	
Pilot Drill, Micro, Short, 9mmL	1.5 Micro; Blue	blue	1.2mmD/9mmL		69.01.14Z
Pilot Drill, Micro, Long, 14mmL	1.5 Micro; Blue	blue	1.2mmD/14mmL	69.01.09Z	
Pilot Block Drill, Micro	1.5 Micro	green	1.6mmD x 15mmL	69.01.16Z	




Fixation Screws 1.5mmD MICRO System	Cartridge (10 screws)	System/Color Indication	Screw Color	Measurements Length mm	Micro 1.5mmD Self drilling	Micro 1.5mmD Self tapping
Screw, self drilling	1	1.5 Micro	green	3.5	68.85.83Z	---
Screw, self drilling	1	1.5 Micro	green	4	68.85.84Z	---
Screw, self drilling	1	1.5 Micro	green	5	68.85.85Z	---
Screw, self drilling	1	1.5 Micro	green	7	68.85.87Z	---
Screw, self tapping	1	1.5 Micro	silver	9	---	68.85.49Z
Screw, self tapping	1	1.5 Micro	silver	11	---	68.85.51Z
Screw, self tapping	1	1.5 Micro	silver	13	---	68.85.53Z
Screw, self tapping	1	1.5 Micro	silver	15	---	68.85.55Z
Screw, self tapping	1	1.5 Micro	silver	17	---	68.85.57Z



2.0mm MINI System	System/Color Indication	Color Marking	Measurements	Kit & Stand Alone	Stand Alone Only
Tray	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	69.01.11Z	
Screw Driver Handle	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	75.23.52Z	
Screw Driver Insert, Short	2.0 Mini; Red	red	17mmL		75.23.21Z
Screw Driver Insert, Long	2.0 Mini; Red	red	55mmL		75.23.22Z
Pilot Drill, Mini, Short	2.0 Mini; Red	red	1.6mmD/9mmL		69.01.15Z
Pilot Block Drill, Mini	2.0 Mini	yellow	2.1mmD x 15mmL		69.01.17Z



Fixation Screws 2.0mm MINI System	Cartridge (10 screws)	System/Color Indication	Screw Color	Measurements Length mm	Micro 2.0mmD Self drilling
Screw, self tapping	1	2.0 Mini; Red	silver	4	68.94.00Z
Screw, self tapping	1	2.0 Mini; Red	silver	5.5	68.94.01Z
Screw, self tapping	1	2.0 Mini; Red	silver	7	68.94.02Z
Screw, self tapping	1	2.0 Mini; Red	silver	9	68.94.03Z
Screw, self tapping	1	2.0 Mini; Red	silver	11	68.94.04Z
Screw, self tapping	1	2.0 Mini; Red	silver	13	68.94.05Z
Screw, self tapping	1	2.0 Mini; Red	silver	15	68.94.06Z
Screw, self tapping	1	2.0 Mini; Red	silver	17	68.94.07Z



Additional Components	Color Indication	Color Marking	Measurements	Stand Alone Only	
Titanium Mesh	Blue & Red	n/a	42x32mm; 0.1mm thick	68.86.95Z	
Titanium Mesh	Blue & Red	n/a	42x32mm; 0.2mm thick	68.86.96Z	
Screw, Rescue	1 Cartridge (10 screws)	Red	pink	7mm	68.94.10Z
Screw, Rescue	1 Cartridge (10 screws)	Red	pink	11mm	68.94.11Z

To learn more about the Zimmer Screw Fixation System, please visit us online at www.zimmerdental.com, or to speak to a sales representative call 1 (800) 854-7019 (in the US) or +1 (760) 929-4300 (outside the US).

This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendation. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.