


Dental Surgical Wash Tray



1 Compartmentalized

Houses small, reusable instruments in two removable compartments, one with grommets to secure the latch-lock end of drills and tools, and a second utility compartment to hold other small instruments.

2 Compatibility

The porous, stainless steel construct allows for detergent and water passage during the cleaning and rinsing processes in automated cleaners and thermo-disinfectors.

3 Universal Application

Holds small instruments for a variety of implant systems provided by Zimmer Biomet Dental.



ZIMMER BIOMET
Your progress. Our promise.®

Dental Surgical Wash Tray

The Dental Surgical Wash Tray is a stainless steel construct which holds small, reusable dental implantology instruments during the automated wash and disinfection process when carried out by a thermo-disinfector.


COMPARTMENTALIZED

- Wash tray contains two removable components
- Grommet compartment may hold up to 50 instruments with ISO-standard latch-lock shank
- Utility compartment with lid holds small, non-latch-lock instruments




COMPATIBILITY

- Stainless steel construct compatible with enzymatic and alkaline solutions used by automated cleaners/thermo-disinfectors
- Porous design allows detergent and water passage during cleaning and rinsing process



UNIVERSAL

- Houses small implantology instruments for a variety of Zimmer Biomet Dental Implant Systems during the automated cleaning process



ORDERING INFORMATION:

Item Number	Description
ZBDWT01	Dental Surgical Wash Tray

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0512 REV A 03/18 ©2018 Zimmer Biomet. All rights reserved.

