


TSV® BellaTek® Encode® Healing Abutment Matrix

D: Platform Diameter
EP: Emergence Profile
H: Height



Grooves on the occlusal surface of the screw identify the implant system. Three grooves indicate the Tapered Screw-Vent® Implant System.

	3.5 mmD		4.5 mmD			5.7 mmD
	3.8 mmEP	5.0 mmEP	5.0 mmEP	5.6 mmEP	6.0 mmEP	6.8 mmEP
3.0 mmH						
	TEHA3383	TEHA3503	TEHA4503	TEHA4563	TEHA4603	TEHA5683
5.0 mmH						
	TEHA3385	TEHA3505	TEHA4505	TEHA4565	TEHA4605	TEHA5685
7.0 mmH		-	-		-	-
	TEHA3387			TEHA4567		

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Tel: +1-561-776-6700
Fax: +1-561-776-1272

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com



Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZBINST0013 REV B 01/21 ©2021 Zimmer Biomet. All rights reserved.

