



SwissPlus[®] Implant System



Product Catalog



zimmer | dental
Confidence in your hands[™]

How to Order

To order, call Zimmer Dental Customer Service between 6:00 a.m. and 4:30 p.m. PST, Monday through Friday, or fax your order to 888 225 2483. Overnight delivery is available, if requested. Overnight delivery orders must be placed by 2:30 p.m. PST.

To make ordering easier, please have your customer number, the catalog number and the quantity of items you wish to order. Prices do not include applicable sales tax, shipping and handling costs. Specifications and prices are subject to change without notification.

Zimmer Dental accepts VISA, Mastercard and American Express.

Payment terms are 2% 10, net 30. Discount pertains to cash payments only. All payments must be in U.S. Dollars (no foreign currency accepted).

To order call:

In the U.S. 800 854 7019

To fax an order 888 225 2483

Outside the U.S. +1 760 929 4300

Return Policy

U.S. customers, please contact Zimmer Dental's Customer Service Department prior to returning product.

Please observe the following guidelines:

1. Customers may return **unopened and undamaged** Zimmer Dental product in its original package only. All returned product must be in saleable condition with the manufacturing seals intact.
2. Returned product must be accompanied by a Customer Return Number obtained from a Customer Service

representative. If product is returned to Zimmer Dental without a Customer Return Number, a 15% restocking fee will be charged.

3. No credits, exchanges or refunds will be provided on discontinued product, product with a shelf life of less than six months, or instruments and tooling over 90 days. Credit or refund for returned items will be issued upon receipt and inspection of goods by Zimmer Dental.
4. All freight must be prepaid on returned merchandise.
5. Customer returning product within **30 days** of purchase may choose from the following options:
 - Full credit to customer account
 - Full exchange on product
 - Full refund
6. Customer returning product within **60 days** of purchase may choose from the following options:
 - Full credit to customer account
 - Full exchange on product
 - Refund minus 15% restocking fee
7. Customers returning product within **61-90** days of purchase may choose from the following options:
 - Credit to customer account minus 15% restocking fee
 - Full exchange on all product
8. After 90 days, only implant bodies can be returned for exchange (up to 120 days after purchase). No other product returns will be accepted after 90 days.

Puros® Allograft Return Policy

Please note: Zimmer Dental does not permit the return of or offer refunds for any of the *Puros* Allograft products.

Return Procedure

Contact Zimmer Dental's Customer Service Department at 800 854 7019 to obtain a Customer Return Number. This number must appear on the outside of the shipping box that contains the items being returned to Zimmer Dental. **Do not mark any of the product packages or labeling.** Any returned package received without a Customer Return Number will be charged a 15% restocking fee. Returns will not be accepted for product that is not in saleable condition. Please send all packages, in a box, to the following address:

Zimmer Dental Inc.
Attn: Returns Department
1900 Aston Avenue
Carlsbad, CA 92008

International customers, please contact your local Zimmer Dental agent or dealer for return policies.

Trademark Information

ImplantMED is a registered trademark of W&H Dentalwerk Bürmoos GmbH. Locator is a registered trademark of and Locator attachments are manufactured by Zest Anchors, Inc. Tutoplast is a registered trademark of Tutogen Medical, Inc. CollaTape, CollaPlug, CollaCote, BioMend, and BioMend Extend are manufactured for Zimmer Dental Inc. by Integra LifeSciences Corporation. CollaTape and CollaPlug are registered trademarks of and CollaCote is a trademark of Integra LifeSciences Corporation.

Table of Contents

Implants

Implant Dimensions	1
Implants	2
Surgical Cover Screws	2

Prosthetics

Components for Impression Transfers	3
Components for Cement-Retained Restorations	4
Components for Custom Restorations	6
Components for Screw-Retained Restorations	7
Overdenture Attachments	9
Prosthetic Tools and Torque Wrenches	12

Surgical Accessories

<i>Zimmer</i> ® Surgical Motor System	13
Surgical Drape Kit	14

Surgical Tools and Information

Surgical Tools	15
Surgical Tray Diagram and Protocol	16-17

Zimmer Regenerative Products

Hard Tissue Products	21
Soft Tissue Products	22

Education Products

Patient Education Material	25
----------------------------	----

Product Index




Index	26
-------	----

Note: Images shown in the catalog may not be to scale.

A Leader in Dental Implant Innovation

Proprietary Technology from Zimmer Dental:

- Internal Hex or Octagon Connection: Prosthetic Simplicity
- Multiple-Lead Threads: Faster Insertion
- Tapered for Bone Expansion: Improves Initial Stability
- Fixture Mount / Transfer Pkg.: Surgical Simplicity / Prosthetic Versatility
- Method for Immediately Placing a Non-Occlusive Dental Implant Prosthesis
- Cutting Drill for Endosseous Implants



Tapered SwissPlus Implants with MTX Surface

Includes Fixture Mount/Abutment & 2mm Long Surgical Cover Screw.



Catalog Numbers						
Implant Diameter	Abutment Platform Diameter	Connection	8mmL	10mmL	12mmL	14mmL
3.7mm	3.8mm	2.5mmD Internal Hexagon	SPMB8	SPMB10	SPMB12	SPMB14
3.7mm	4.8mm	3.0mmD Internal Octagon	SPB8	SPB10	SPB12	SPB14
4.8mm	4.8mm	3.0mmD Internal Octagon	SPWB8	SPWB10	SPWB12	SPWB14

Straight SwissPlus Implants with MTX Surface

Includes Fixture Mount/Abutment & 2mm Long Surgical Cover Screw.



Catalog Numbers						
Implant Diameter	Abutment Platform Diameter	Connection	8mmL	10mmL	12mmL	14mmL
4.1mm	4.8mm	3.0mmD Internal Octagon	OPB8	OPB10	OPB12	OPB14
4.8mm	4.8mm	3.0mmD Internal Octagon	OPWB8	OPWB10	OPWB12	OPWB14

Surgical Cover Screws

Catalog Numbers					
	Platform Diameter	Component Length (L)	Component Diameter (ø)	Catalog No.	
	Narrow Diameter Implant Short Surgical Cover Screw	3.8mm	1mm	4.6mm	SPSCS
	Narrow Diameter Implant Surgical Cover Screw	3.8mm	2mm	4.6mm	SPSC
	Narrow Diameter Implant Long Surgical Cover Screw	3.8mm	3mm	4.6mm	SPSCL
	Short Surgical Cover Screw	4.8mm	1mm	5.2mm	OPSCS
	Surgical Cover Screw	4.8mm	2mm	5.2mm	OPSC
	Long Surgical Cover Screw Includes Suture Groove	4.8mm	3mm	5.2mm	OPSCL
	Wide Diameter Surgical Cover Screw	4.8mm	2mm	6.0mm	OPWSC
	Long Wide Diameter Surgical Cover Screw Includes Suture Groove	4.8mm	3mm	6.0mm	OPWSCL

Easy Impressions


The *SwissPlus* Plastic Transfer is used with the Straight Abutment for a press-on impression of the abutment. Align and press-fit the Castable Plastic Coping to the abutment. The retentive plug contacts the flat of the abutment for rotational stability.



Continue with your impression. The matching abutment analog duplicates the shoulder of the implant and abutment. The coping can also be used to create a reduction coping if the abutment is prepped.

The Straight Abutment and Plastic Transfer are available for both the 3.8mmD and 4.8mmD Implant Platforms.


Impression-Taking Flow Chart



Transfer Components

Catalog Numbers				
	Description	Implant Platform Diameter	Emergence Profile Diameter	Catalog No.
	Indirect (Closed-Tray) Transfer – Includes (OPTS)	4.8mm	5.2mm	OPT/5
	Indirect (Closed-Tray) Transfer – Includes (OPTS)	4.8mm	6.0mm	OPT/6
	Indirect (Closed-Tray) Transfer for Internal Hexagon – Includes (HLTSI)	3.8mm	4.6mm	SPMT
	Implant Replica with Internal Octagon	4.8mm	•	OPR
	Implant Replica with Internal Hexagon	3.8mm	•	SPMR
	Narrow Abutment/Shoulder Replica	4.8mm	•	OPAR
	Narrow Abutment/Shoulder Replica	3.8mm	•	SPMAR

Cement-Retained Restoration Flow Chart



Straight Abutments

Catalog Numbers				
	Description	Implant Platform Diameter	Emergence Profile Diameter	Catalog No.
	Narrow Straight Abutment (One-Piece) Includes Plastic Transfer/Waxing Coping (OPC), Does Not Engage Internal Octagon	4.8mm	3.5mm	OPA
	Narrow Straight Abutment (Two-Piece) Includes Plastic Transfer/Waxing Coping (OPC) and Screw (GPCAS), Engages Internal Octagon	4.8mm	3.5mm	OPA2
	Narrow Straight Abutment (Two-Piece) Includes Plastic Transfer/Waxing Coping (SPMC) and Screw (MHLAS), Engages Internal Hexagon	3.8mm	2.9mm	SPMA2
	Replacement Plastic Waxing Coping	4.8mm	•	OPC
	Replacement Plastic Waxing Coping	3.8mm	•	SPMC
	Straight Abutment/Transfer Includes Screw (OPAS), Engages Internal Octagon	4.8mm	5.2mm	OPA/5
	Straight Abutment/Transfer Includes Screw (OPAS), Engages Internal Octagon	4.8mm	6.0mm	OPA/6
	Straight Abutment/Transfer Includes Screw (URS3), Engages Internal Hexagon	3.8mm	4.6mm	FMTM2
	Replacement Screw for Narrow Straight Abutment (OPA2)	4.8mm	•	GPCAS
	Replacement Screw for Narrow Straight Abutment (SPMA2)	3.8mm	•	MHLAS
	Replacement Transfer Screw for Straight Abutment/Transfer (OPA/5 and OPA/6)	4.8mm	•	OPAS
	Replacement Transfer Screw for Straight Abutment/Transfer (FMTM2)*	3.8mm	•	URS3

* Call for availability.

Multi-Functional Fixture Mount/Transfer

SwissPlus Implants are value-packaged with a Fixture Mount that functions as an impression post and preparable final abutment.



Implant Insertion

Fixture Mount Drill used to insert self-tapping implant.



Transfer Impression


Block out the top of the Fixture Mount/Transfer to make a conventional impression.



Preparable Final Abutment

Shorten and prepare the Fixture Mount/Transfer for use as an abutment.

Cement - Retained Angled Restoration Flow Chart



Straight Abutments


Catalog Numbers				
	Description	Implant Platform Diameter	Emergence Profile Diameter	Catalog No.
	Fixture Mount/Transfer Regular Includes Screw (URS1)	4.8mm	5.2mm	FMS
	Fixture Mount/Transfer Wide Includes Screw (URS1)	4.8mm	6.0mm	FMSW
	Fixture Mount/Transfer Mini Includes Screw (URS2)	3.8mm	4.6mm	FMSM
	Fixture Mount Screw (Regular and Wide)*	•	•	URS1
	Fixture Mount Screw (Mini)*	•	•	URS2

* Call for availability.

20° Angled Abutments

Catalog Numbers				
	Description	Implant Platform Diameter	Emergence Profile Diameter	Catalog No.
	20° Angled Abutment for Internal Octagon Includes Screw (OPH20S)	4.8mm	5.5mm	OPH20
	Replacement Screw for 20° Angled Abutment	•	•	OPH20S
	20° Angled Abutment for Internal Hexagon Includes Screw (AH20S)	3.8mm	4.6mm	SPH20
	Replacement Screw for 20° Angled Abutment	•	•	AH20S

Custom Restoration Flow Chart



“Cast-To” Gold Abutments, Engaging

Catalog Numbers				
	Description	Implant Platform Diameter	Emergence Profile Diameter	Catalog No.
	Gold Abutment Includes Plastic Sheath & Screw (OPS & GPCAS)	4.8mm	5.2mm	OPGA
	Replacement Screw for “Cast-To” Abutments, Copings & Castable Copings	4.8mm	•	GPCAS
	Waxing Screw for “Cast-To” Abutments, Copings & Castable Copings	4.8mm	•	WSX
	Gold Abutment Includes Plastic Sheath & Screw (OPS & MLHAS)	3.8mm	4.6mm	SPGA
	Replacement Screw for Narrow “Cast-To” Abutments	3.8mm	•	MHLAS
	Waxing Screw for “Cast-To” Abutments	3.8mm	•	MTWSD
	Press-Fit Castable Plastic Sheath for “Cast-To” Abutments and Copings	•	•	OPS

“Cast-To” Gold & Plastic Copings, Non-Engaging

Used for implant-level, multiple-unit, screw-retained restorations. Includes Plastic Sheath & Screw (OPS & GPCAS).

Catalog Numbers				
	Description	Implant Platform Diameter	Emergence Profile Diameter	Catalog No.
	Gold Coping Includes Plastic Sheath & Screw (OPS & GPCAS)	4.8mm	5.0mm	OPGC
	Plastic Castable Coping Includes Screw (GPCAS)	4.8mm	5.0mm	OPCC


The Zimmer Institute™



The Zimmer Institute (ZI) Dental Programs are designed and dedicated to helping practitioners obtain the knowledge, skills, and confidence that are essential for the practice of contemporary implant dentistry. The Zimmer Institute Programs – either continuing education or Zimmer product training – are developed and presented by outstanding university faculty and private practitioners from around the world and are currently divided into four levels of instruction.

These programs provide both general practitioners and specialists with outstanding educational experiences for learning today's implant technology and, in the future, will provide an avenue for learning and incorporating "tomorrow's implant technology" as well. In essence, they help to ensure the very best in patient care by providing practitioners with an exceptional pathway for lifetime learning.


Screw - Retained Restoration Flow Chart



Tapered Abutments

Includes Healing Cap (TATHC).

Catalog Numbers

Description	Implant Platform Diameter	Cuff Height	Catalog No.
 Tapered Abutment	4.8mm	1.6mm	OPACT
Tapered Abutment	4.8mm	3.0mm	OPACT3

Tapered Healing Cap

Healing Cap threads onto the Tapered Abutment.




Catalog Numbers

Description	Catalog No.
 Titanium Healing Cap	TATHC

Tapered Abutment Transfer Components

Catalog Numbers		
	Description	Catalog No.
	Direct Transfer (Open-Tray) Includes Transfer Screw (SCDTS)	ACTDT
	Indirect Transfer (Closed-Tray)	ACTIT
	Abutment Replica	ACTR

Tapered Abutment Copings, Sheaths & Bar Systems

Catalog Numbers			
	Description	Length	Catalog No.
	Tapered Abutment Gold Coping – Includes Sheath (OPS) and Screw (SCTS)	•	ACTGC
	Plastic Castable Coping Includes Screw (SCTS)	5.0mm	ACTP
	Bar Gold Coping – Includes Screw (SCTS)	3.0mm	TGC3
	Bar Gold Coping – Includes Screw (SCTS)	5.0mm	TGC5
	Titanium Bar Coping Includes Screw (SCTS)	5.0mm	TTC5
	Titanium Temporary Coping Includes Screw (SCTS)	9.5mm	ACTT
	Replacement Screw for Copings	•	SCTS
	Abutment Waxing Screw	•	SCWS
	1.8mm Round Gold Bar Includes 3 <i>Hader</i> [®] Bar Clips, 3 Processing Clips, 1 Insertion Tool	50mm	HGB
	Dolder Resilient Gold Bar Includes Machined Metal Sleeve and Spacer	50mm	DGB
	<i>Hader</i> Bar Clip System Includes 6 Processing Clips, 6 <i>Hader</i> Bar Clips, 6 <i>Hader</i> Bar Patterns, 1 Insertion Tool	•	BS1

Locator® Overdenture Attachments

Locator Attachments are available for the SwissPlus Dental Implant System.



The *Locator* Overdenture Attachment is the ideal solution for overdenture restorations.


- **Low Profile.** A low vertical height makes the *Locator* Overdenture Attachment an effective option for treatment of a wider range of patients.
- **Easy for Patients.** Self-locating attachment gives patients the ability to easily position and place their overdentures.
- **Predictable.** Dual-retention and tremendous durability ensure long-lasting stability and performance.¹
- **Versatile.** A variety of options for retention and angulation correction lets you tailor the case to fit the needs of the patient.

The *Locator* Bar Attachment System is available for use when the splinting of implants is desired to increase the stability of a bar. This system encompasses the same self-aligning feature, dual-retention, long-lasting durability, and low attachment profile as the *Locator* Overdenture Attachment System.

The *Locator* Bar Attachment System provides three options for the fabrication of a resilient attachment on an implant-retained cast alloy or milled titanium bar: a Drill and Tap system, a Laser Welding system, and a Cast-To option.

¹Data on file.

Locator Overdenture Restoration Flowchart





Locator Abutments for SwissPlus Implants

Catalog Numbers				
		Platform Diameter	Cuff Height	Catalog No.
	Locator Abutment	4.8mm	1.29mm	OPLOC4/1
	Locator Abutment	4.8mm	2mm	OPLOC4/2
	Locator Abutment	4.8mm	3mm	OPLOC4/3
	Locator Abutment	4.8mm	4mm	OPLOC4/4
	Locator Abutment	3.8mm	0.75mm	OPLOC3/1
	Locator Abutment	3.8mm	2mm	OPLOC3/2
	Locator Abutment	3.8mm	3mm	OPLOC3/3
	Locator Abutment	3.8mm	4mm	OPLOC3/4


Locator Attachment Components

Catalog Numbers			
	Description	Qty	Catalog No.
	Male Processing Package, Includes 2 Each of Metal Housings with Black Processing Males, Block-Out Spacers, Clear, Pink and Blue Replacement Males	2 Ea	LOCMP2
	Male Processing Package, Includes 10 Each of Parts Listed Above	10 Ea	LOCMP10
	Paralleling Post	4	LOCPP
	Replacement Male, 5 lbs. (Clear)	4	LOCRMW
	Replacement Male, Light Retention, 3 lbs. (Pink)	4	LOCRMP
	Replacement Male, Extra Light Retention, 1.5 lbs. (Blue)	4	LOCRMB
	Extended Range Replacement Male, 4 lbs. (Green)	4	LOCRMG
	Extended Range Replacement Male, Light Retention, 1.5 lbs. (Red)	4	LOCRMW


Locator Abutment Transfer Components

Catalog Numbers			
	Description	Qty	Catalog No.
	Impression Coping, Includes Black Processing Male	4	LOCIMP
	Female Analog, 4mmD (For Use with 3.8mmD <i>SwissPlus</i> Implant)	4	LOCAN/4
	Female Analog, 5mmD (For Use with 4.8mmD <i>SwissPlus</i> Implant)	4	LOCAN/5


Locator Bar Attachment System – Drill and Tap Components

Catalog Numbers			
	Description	Qty	Catalog No.
	<i>Locator</i> Bar Female (2mm Thread)	2	LOCBF2
	<i>Locator</i> Yellow Bar Processing Package	2	LOCBPP
	<i>Locator</i> Bar Drill and Tap Holder	1	LOCPTH
	<i>Locator</i> 2.0mm Castable Threaded Insert (Optional)	10	LOCCTI
	<i>Locator</i> 2.0mm Bar Tap (2mm Thread)	1	LOCBT2
	<i>Locator</i> 1.7mm Bar Drill (2mm Thread)	1	LOCBD





Locator Bar Attachment System – Laser Weld Components

Catalog Numbers			
	Description	Qty	Catalog No.
	<i>Locator</i> Laser Bar Female (Stainless Steel)	2	LOCLBF
	<i>Locator</i> Yellow Bar Processing Package	2	LOCBPP
	<i>Locator</i> Paralleling Mandrel	1	LOCPMAND

Locator Bar Attachment System – Cast-To Components

Catalog Numbers			
	Description	Qty	Catalog No.
	<i>Locator</i> Cast-To Bar Female	2	LOCCTB
	<i>Locator</i> Yellow Bar Processing Package	2	LOCBPP
	<i>Locator</i> Paralleling Mandrel	1	LOCPMAND

Locator Tools

Catalog Numbers			
	Description	Qty	Catalog No.
	<i>Locator</i> Core Tool	1	LOCCT2
	<i>Locator</i> Core Tool Replacement Tip (Fits LOCCT)	1	LOCCTR
	<i>Locator</i> Torque Wrench Insert Driver, 15mmL	1	LOCTW15
	<i>Locator</i> Torque Wrench Insert Driver, 21mmL	1	LOCTW21

Platforms - Internal Octagon and Internal Hexagon

Two Internal Platforms


Most *SwissPlus* Implants feature a 4.8mmD platform with an internal octagon. A full range of prosthetic components is available for simple esthetic restorations.



The 3.8mm platform features an internal hexagon and is only available on the 3.7mm *Tapered SwissPlus* Implant. The prosthetic options for this implant include cementable and custom abutments, which are well-suited for anterior esthetic restorations.



Ball and Cap Overdenture Restoration Flow Chart



Ball Abutments

Includes Cap Attachment Housing and Nylon Liner.

Catalog Numbers

	Description	Platform Diameter	Cuff Height	Catalog No.
	Ball Abutment & Cap Attachment	4.8mm	1.6mm	OPBA
	Ball Abutment & Cap Attachment	4.8mm	3.0mm	OPBA3

Ball Abutment Transfers & Replicas

Catalog Numbers

	Description	Catalog No.
	Ball Abutment Transfer (Package of 2)	OPBAT
	Ball Abutment Replica	OPBAR

Overdenture Attachment Prosthetic Components

Catalog Numbers

	Description	Catalog No.
	Cap Attachment System Includes 4 Nylon Liners, 2 Positioning Rings, 4 Housings, 4 Castable Ball Patterns	CAS
	Cap Attachment Instruments Includes Seating Tool, Reamer, Paralleling Tool	CAI
	Cap Attachment Housing & Cap Attachment Nylon Liner (CAH & CAN)	CA
	Cap Attachment Housing	CAH
	Cap Attachment Nylon Liner (Transparent)	CAN
	Cap Attachment Nylon Liner (Grey — Rigid Retention)	CAN-G
	Cap Attachment Transfer (Yellow)	CAT
	Castable Ball Pattern (2 Balls per Pattern)	CAB
	Micro Cap Attachment Nylon Liner	CANM

Prosthetic Tools and Torque Wrenches

Catalog Numbers		
	Description	Catalog No.
	Hex Tool for Surgical and Prosthetic Screws (1.25mmD, 17mmL)	HX1.25
	Long Hex Tool for Surgical and Prosthetic Screws (1.25mmD, 22mmL)	HXL1.25
	Stainless Steel Tapered Frictional Hex Tool (1.25mmD, 23mmL)	THX1.25
	Long Stainless Steel Tapered Frictional Hex Tool (1.25mmD, 27mmL)	THXL1.25
	Hex Drill for Surgical and Prosthetic Screws (Requires Reduction Handpiece)	HX1.25D
	Torque Wrench, 30 Ncm	TW30
	Torque Wrench, 20 Ncm	TW20
	Torque Wrench Hex Driver 1.25mmD, for TW20 and TW30	TW1.25
	Long Torque Wrench Hex Driver 1.25mmD, for TW20 and TW30	TW1.25L
	Locator Torque Wrench Insert Driver, 15mmL, for TW20 and TW30	LOCTW15
	Locator Torque Wrench Driver, 21mmL, for TW20 and TW30	LOCTW21
	Locator Core Tool	LOCCT2
	Locator Core Tool Replacement Tip (Fits LOCCT)	LOCCTR
	Prosthetic Kit Tray, Autoclavable (Empty)	2320
	PTS Disposable Inserts, 20 Ncm, White (Package of 5)	PTS20
	PTS Disposable Inserts, 30 Ncm, Black (Package of 5)	PTS30
	Abutment Holder	ABTH
	Reamer for Tapered Abutment Copings and <i>SwissPlus</i> Copings and Abutments	PR
	Removal Tool for Abutment Screws or Fixation Screws with 1.25mmD Hex, Carbide Steel	SRT
	Removal Tool for Internal Hex Implant Interface	TLRT2
	Removal Tool for Internal Octagon Implant Interface	HLRTX2

Dríva™ Surgical Drills

High performance *Dríva* Reusable Drills offer clear-cut advantages.

Through extensive research and testing, Zimmer Dental has developed the *Dríva* Reusable Drill, providing high levels of strength and cutting efficiency. The progressive design improves confidence when creating an osteotomy by delivering:

- Improved coating. Revolutionary corrosion-resistant coating protects drill sharpness.
- Enhanced visibility. Grooved and laser-etched depth bands facilitate visual measurement, along with a matte finish that reduces glare, for an optimal surgical experience.



Zimmer Surgical Motor System

System Features

- Designed for implantology and other maxillofacial applications
- Simple, user-friendly programs and maintenance-free features for convenience and time savings
- Increased torque control for confidence during surgeries
- Quiet performance of motor and handpieces
- Wide speed range of 300 to 40,000 rpm (1:1)
- Worldwide service locations**



Zimmer Surgical Motor System

Catalog Numbers

Description	Catalog No.	Qty
Zimmer Surgical Motor System* (100-130V)	00900125	1
Zimmer Surgical Motor System* (220-240V)	00900124	1

* Each system includes: console, motor, motor holder, foot control, irrigation stand, testing gauge for burs, oil spray cap for handpieces, 3 sets of one-piece, disposable irrigation tubing and power cords.

** Contact your local sales representative, distributor or call Zimmer Dental Customer Service for information regarding service center locations.

Surgical Handpieces

Catalog Numbers

Description	Catalog No.	Qty
WI-75 E/KM Contra-Angle Handpiece 20:1 (Includes Y-Connector Assembly, Internal Bur Coolant Tube)	10207557	1
WS-92 Contra-Angle Handpiece 1:2	10209201	1
S-11 Straight Handpiece 1:1	00001100	1
S-12 Angled Handpiece 1:2	10101200	1
S-8 R Reciprocating Saw	10100801	1
S-8 O Oscillating Saw	10100802	1
S-8 S Sagittal Saw	10100803	1

Surgical Saw Blades

Catalog Numbers

Description	Length	Catalog No. (1/Pkg)	Catalog No. (5/Pkg)
Reciprocating	15mm	04726000	04800100
Reciprocating	20mm	04725900	04800000
Reciprocating	30mm	04725700	04799900
Oscillating	8mm	04800500	N/A
Oscillating	15mm	04800400	N/A
Sagittal	6mm	04726200	04800300
Sagittal	10mm	04726100	04800200

ATTENTION: Part #02139800, Testing gauge for burs provided with the motor is a required part for proper use of all W&H handpieces. This gauge tests the condition of the burs for wear. If burs are not tested prior to use, they can jam inside the handpiece and severely damage the mechanisms. Call Zimmer Dental for additional information.

Irrigation Accessories

Catalog Numbers		
Description	Catalog No.	Qty
One-Piece Irrigation Tubing Set, Disposable	04363600	6
One-Piece Irrigation Tubing Set, Autoclave	04719400	1
Irrigation Spike with Roller, Autoclave	0285540R	1
Pump Tubing with Connectors, Autoclave	04013900	1
Spare Pump Tubing without Connectors, Autoclave	04014000	3
Irrigation Tubing, 2.2m, Disposable	00929300	10
Internal Bur Coolant Tube	02610500	1
Y-Connector	02653901	1
Y-Tube Assembly	02654700	1
Set of Tube Clamps	04019000	1
Spray Clip for External/Internal Irrigation	04757100	1
Irrigation Stand	04005900	1

Micromotors and Accessories

Catalog Numbers		
Description	Catalog No.	Qty
Motor with Cable for Surgical Motor System, 1.8m	04720022	1
Motor with Cable for <i>ImplantMED</i> ® Surgical Motors, 1.8m	04009600	1
Motor Holder for <i>Zimmer</i> Surgical Motor System	04735700	1
Oil Spray Cap for Handpieces*	02038200	1
Testing Gauge for Burs	02139800	1
Sterilizing Protector for Motor	04032600	1
Sterilizing Box for Motor	04013500	1
Transportation Case	04013600	1
ST-GC Treatment Trolley with Shelf	04541900	1

* Call for service oil ordering instructions.

Aseptic Surgical Drape Kit

Catalog Numbers		
Model	Catalog No.	Qty
Complete Surgical Drape Kit	ASD	1 Ea
Split Drape (53"x 77")	Hose Cover	
Towels (2)	Tray Foam (14"x 18")	
Light Handle Covers (2)	Table Cover	
Sterile Pouch	CSR Wrap	
Cap	Foam (2)	
Gowns (2)		

Tapered SwissPlus Implant



1. Fixture Mount

is both a transfer for fixture-level impressions and a full-contour preparable abutment that can allow margins to follow tissue contour.

2. Tapered SwissPlus Implant

features double-lead threads for faster insertion and a tapered body with an osteocompressive surgical protocol designed to enhance initial stability in soft bone!

3. MTX Surface

increases surface texture for enhanced bone apposition?

4. Self-Tapping Design

means less surgical time.


^{1,2} Data on file.


SwissPlus Implant Surgical Kit

Catalog Numbers

Description	Catalog No.	Qty/Pkg
Complete Surgical Kit	OPCST	1 Ea
Surgical Tray	OPSUR	
2.3mmD <i>Dríva</i> Pilot Drill, 17mmL	OP2.3D	
2.3mmD <i>Dríva</i> Pilot Drill, 11mmL	SV2.3DS	
2.8mmD <i>Dríva</i> Intermediate Drill, 17mmL	OP2.8D	
2.8mmD <i>Dríva</i> Intermediate Drill, 11mmL	SV2.8DS	
3.4mm/2.8mmD <i>Dríva</i> Final Step Drill, 17mmL	SP3D	
3.4mm/2.8mmD <i>Dríva</i> Final Step Drill, 11mmL	TSV3DS	
3.5mmD <i>Dríva</i> Drill, 17mmL	OP3.5D	
3.5mmD <i>Dríva</i> Drill, 11mmL	TLSDS	
4.4mm/3.8mmD <i>Dríva</i> Final Step Drill, 17mmL	SP4D	
4.4mm/3.8mmD <i>Dríva</i> Final Step Drill, 11mmL	TSV4DS	
4.2mmD <i>Dríva</i> Drill, 17mmL	OP4.2D	
4.2mmD <i>Dríva</i> Drill, 11mmL	SVWDS	
Removal Tool/Long Tapered 1.25mmD Hex Tool	HLRTX2	
Tapered 1.25mmD Hex Tool	THX1.25	
1.25mmD Hex Tool, 22mmL	HXL1.25	
1.25mmD Hex Tool, 17mmL	HX1.25	
1.25mmD Hex Drill	HX1.25D	
2.5mmD <i>GemLock</i> ® Retaining Hex Driver, Short, 17mmL	RH2.5	
2.5mmD <i>GemLock</i> Retaining Hex Driver, Long, 28mmL	RHL2.5	
2.5mmD <i>GemLock</i> Retaining Hex Drill, 21mmL	RHD2.5	
3.0mmD Octagon Insertion Tool, Short, 17mmL	OT3.0-S	
3.0mmD Octagon Insertion Tool, Long, 29mmL	OTL3.0-S	
Drill Extender	DE	
Paralleling Tool (4 Pieces)	PPAR	
<i>GemLock</i> Retaining Square Ratchet	RSR	
Stainless Steel Screwdriver Handle with Square Connection	SSHS	
Surgical Tray	OPSUR	1 Ea
Tray Only – Does Not Include Instruments		

Note: For detailed image of tray see page 16.






Drilling Sequence for Tapered SwissPlus Implant


3.7mmD Tapered SwissPlus Implant (3.8mmD Platform and 4.8mmD Platform)



3.7mmD



1
OP2.3D/SV2.3DS
 2.3mmD
 Drill



2
FOR SOFT BONE
OP2.8D/SV2.8DS
 2.8mmD
 Drill



2
FOR DENSE BONE
SP3D/TSV3DS
 3.4mm/2.8mmD
 Drill



3
OPTIONAL FOR
DENSE BONE
DT3.7
 Cortical Tap

4.8mmD Tapered SwissPlus Implant (4.8mmD Platform)



4.8mmD



1
OP2.3D/SV2.3DS
 2.3mmD
 Drill



2
SP3D/TSV3DS
 3.4mm/2.8mmD
 Drill



3
FOR SOFT BONE
OP3.5D/TLSDS
 3.5mmD
 Drill



3
FOR DENSE BONE
SP4D/TSV4DS
 4.4mm/3.8mmD
 Drill



4
OPTIONAL FOR
DENSE BONE
DT4.8
 Cortical Tap

Drilling Sequence for Straight SwissPlus Implant

4.1mmD Straight SwissPlus Implant (4.8mmD Platform)



4.1mmD



1
OP2.3D/SV2.3DS
 2.3mmD
 Drill



2
OP2.8D/SV2.8DS
 2.8mmD
 Drill



3
OP3.5D/TLSDS
 3.5mmD
 Drill

4.8mmD Straight SwissPlus Implant (4.8mmD Platform)



4.8mmD



1
OP2.3D/SV2.3DS
 2.3mmD
 Drill



2
OP2.8D/SV2.8DS
 2.8mmD
 Drill







3
OP3.5D/TLSDS
 3.5mmD
 Drill




4
OP4.2D/SVWDS
 4.2mmD
 Drill



Dríva Drills

Catalog Numbers			
Diameter		11mmL	17mmL
 2.3mm		SV2.3DS	OP2.3D
 2.8mm		SV2.8DS	OP2.8D
 3.5mm		TLSDS	OP3.5D
 4.2mm		SVWDS	OP4.2D


Optional Drill

Catalog Numbers		
Diameter		8mmL
 Drill 2.1mm/1.6mm, 8mmL		0201



Dríva Step Drills

Catalog Numbers			
Diameter		11mmL	17mmL
 3.4mm/2.8mm		TSV3DS	SP3D
 4.4mm/3.8mm		TSV4DS	SP4D

Cleaning Tool

Catalog Numbers			
Description		Catalog No.	Qty/Pkg
 Dríva Drill Cleaning Wire		NM1940	1 Ea


Cortical Bone Tap for Tapered SwissPlus Implants

Catalog Numbers		
Description		Catalog No.
 3.7mm Double-Lead Cortical Tap Tool		DT3.7
 4.8mm Double-Lead Cortical Tap Tool		DT4.8

MTX Surface Allows for Increased Bone Apposition¹

Zimmer's proprietary *MTX* Surface has been shown to allow for increased bone apposition compared to machined surfaces¹. In a recent clinical evaluation the *MTX* Surface achieved excellent bone-to-implant contact and osteoconductive capacity².

- Implants are textured with soluble, biocompatible grit-blasting medium
- Implant threads are designed for maximum cutting efficiency
- Cutting grooves remain intact for efficient self-tapping
- *MTX* process results in a uniform texture of closely-spaced micropits








¹Trisi P, Marcato C, Todisco M. Bone-to-implant apposition with machined and MTX micro-textured implant surfaces in human sinus grafts. *Int J Periodontics Restorative Dent.* 2003;23(5):427-437.

²Todisco M, Trisi P. Histomorphometric evaluation of six dental implant surfaces after early loading in augmented human sinuses. *J Oral Implantol* 2006;32(4):153-166.






Single Patient Drills for Tapered SwissPlus

Catalog Numbers		
Description	Contents	Catalog No.
Drill Kit, Single Patient 3.7mm <i>Tapered SwissPlus</i> Implant	3.0mm Round Bur 2.3mm Pilot Drill, 17mmL 3.4mm/2.8mm Stepped Drill, 17mmL	DSSK
Drill Kit, Single Patient 3.7mm <i>Tapered SwissPlus</i> Implant	3.0mm Round Bur 2.3mm Pilot Drill, 11mmL 3.4mm/2.8mm Stepped Drill, 11mmL	DSTK3S
Drill, Single Patient 4.8mm <i>Tapered SwissPlus</i> Implant	4.4mm/3.8mm Stepped Drill, 17mmL	DSSW
Drill, Single Patient 4.8mm <i>Tapered SwissPlus</i> Implant	4.4mm/3.8mm Stepped Drill, 11mmL	DST4S












Surgical Sequence Requirements

3.7mm Tapered Implant				
		DSSK 8, 10, 12 or 14mmL		
		DSTK3S 8 or 10mmL		
4.8mm Tapered Implant				
		DSSK 8, 10, 12 or 14mmL		
		DSTK3S 8 or 10mmL		
		DST4S 8 or 10mmL		


Hex Tools

Catalog Numbers		
	Description	Catalog No.
	Hex Tool for Surgical and Prosthetic Screws (1.25mmD, 17mmL)	HX1.25
	Long Hex Tool for Surgical and Prosthetic Screws (1.25mmD, 22mmL)	HXL1.25
	Tapered Frictional Hex Tool (1.25mmD, 23mmL)	THX1.25
	Tapered Frictional Hex Tool, Long (1.25mmD, 27mmL)	THXL1.25
	Hex Drill for Surgical and Prosthetic Screws (1.25mmD) (Requires Reduction Handpiece)	HX1.25D

Surgical Armamentaria

Catalog Numbers		
	Description	Catalog No.
	Screwdriver Handle	SSHS
	<i>GemLock</i> Retaining Square Ratchet	RSR
	Drill Extender	DE
	Paralleling Tool	PPAR
	Octagon Insertion Tool For Implants with 3.0mmD Octagon (17mmL)	OT3.0-S
	Octagon Insertion Tool, Long For Implants with 3.0mmD Octagon (29mmL)	OTL3.0-S
	<i>GemLock</i> Retaining Hex Drill (2.5mmD, 21mmL)	RHD2.5
	<i>GemLock</i> Retaining Hex Driver, Short (2.5mmD, 17mmL)	RH2.5
	<i>GemLock</i> Retaining Hex Driver, Long (2.5mmD, 28mmL)	RHL2.5
	Removal Tool for Screw-Type Implants, Carbide Steel (Use Cold Sterilization Only – Do Not Autoclave)	IRT
	Removal Tool for Internal Octagon Implant Interface	HLRTX2

Working together with nature.



Filling Nature's Void.

Zimmer Dental offers a complete line of complementary regenerative products to fulfill your clinical needs. Using collagen-based materials derived from natural sources, these products work with the host tissue to regenerate hard and soft tissue structures.

Hard Tissue Products

Zimmer's *Puros* family of allograft bone products offers the clinician excellent clinical results in cases such as sinus elevation, ridge augmentation, ridge preservation, socket repair and periodontal defects.

Soft Tissue Products

Zimmer's family of soft tissue regeneration products includes *Puros* Dermis Allograft Tissue Matrix, *Zimmer* Socket Repair Membrane, *Puros* Pericardium Allograft Barrier Membranes, *BioMend*® Bovine Collagen Barrier Membranes, and Colla Bovine Collagen Wound Dressings. These products are used in conjunction with *Puros* Bone Grafts for Guided Bone and Guided Tissue Regeneration, as well as soft tissue augmentation and healing.

Puros Allografts Refined Using the Unique Tutoplast® Process

The *Puros* Allograft family of products is refined using the proprietary, 5-stage *Tutoplast* Process. The *Tutoplast* Process has been fully validated for preservation of tissue biomechanical properties, viral inactivation, prion inactivation, healing and elimination of antigenicity.



1. Delipidization



2. Osmotic treatment



3. Oxidative treatment



4. Solvent dehydration



5. Limited-dose gamma irradiation

There is one slight difference between the hard and soft tissue processing. *Hard* tissue processing utilizes a delipidization step to inactivate lipid-enveloped viruses and to remove fat, enhancing the healing characteristics. *Soft* tissue processing mechanically removes fat and includes an *alkaline* treatment to reduce prion infectivity.

The *Tutoplast* Process preserves the natural collagen architecture of the graft tissue, providing an excellent matrix for host tissue regeneration and assimilation.



Hard Tissue Products

Puros Bone Allografts


Catalog Numbers		
Description	Catalog No.	
	<i>Puros</i> Cancellous Particles, 0.5cc, 250-1000	8210R
	<i>Puros</i> Cancellous Particles, 1.0cc, 250-1000	8211R
	<i>Puros</i> Cancellous Particles, 2.0cc, 250-1000	8209R
	<i>Puros</i> Cancellous Particles, 0.5cc, 1000-2000	8212R
	<i>Puros</i> Cancellous Particles, 1.0cc, 1000-2000	8213R
	<i>Puros</i> Cancellous Particles, 2.0cc, 1000-2000	8214R
	<i>Puros</i> Cortical Particles, 0.5cc, 250-1000	8271R
	<i>Puros</i> Cortical Particles, 1.0cc, 250-1000	8272R
	<i>Puros</i> Cortical Particles, 2.0cc, 250-1000	8273R
	<i>Puros</i> Cortical Particles, 0.5cc, 1000-2000	8274R
	<i>Puros</i> Cortical Particles, 1.0cc, 1000-2000	8275R
	<i>Puros</i> Cortical Particles, 2.0cc, 1000-2000	8276R
	<i>Puros</i> Block Allograft, 10mm	8220R
	<i>Puros</i> Block Allograft, 15mm	8221R

Soft Tissue Products


Puros Membrane Allografts

Catalog Numbers		
Description	Catalog No.	
	<i>Puros Pericardium Membrane, 15x20mm</i>	8770R
	<i>Puros Pericardium Membrane, 20x30mm</i>	8771R
	<i>Puros Pericardium Membrane, 30x40mm</i>	8772R
	<i>Puros Dermis Tissue Matrix, 10x10mm</i>	8793R
	<i>Puros Dermis Tissue Matrix, 10x20mm</i>	8790R
	<i>Puros Dermis Tissue Matrix, 10x40mm</i>	8791R
	<i>Puros Dermis Tissue Matrix, 20x40mm</i>	8792R

Zimmer Socket Repair Membrane

Catalog Numbers		
Description	Catalog No.	
	<i>Zimmer Socket Repair Membrane, 10x20mm</i>	0154

BioMend and BioMend® Extend™ Absorbable Collagen Membranes




Catalog Numbers		
Description	Catalog No.	
	<i>BioMend Membrane, 15x20mm</i>	0103
	<i>BioMend Membrane, 20x30mm</i>	0105
	<i>BioMend Membrane, 30x40mm</i>	0107
	<i>BioMend Extend Membrane, 15x20mm</i>	0140
	<i>BioMend Extend Membrane, 20x30mm</i>	0141
	<i>BioMend Extend Membrane, 30x40mm</i>	0142

The Collagen Advantage



Colla Wound Dressings, *BioMend* Membranes, and *Zimmer Socket* Repair Membrane are composed of Type I Collagen derived from bovine tendon. Due to the physiochemical properties of the collagen matrix, these products promote hemostasis and wound healing, are nonimmunogenic, and are resorbed into the body.

CollaTape®, CollaCote™ & CollaPlug® Collagen Wound Dressings

Catalog Numbers	
Description	Catalog No.
 CollaTape Wound Dressing, 10/Pkg	0100
 CollaCote Wound Dressing, 10/Pkg	0101
 CollaPlug Wound Dressing, 10/Pkg	0102

Calcitite® HA 2040 Particles – Nonresorbable Dense Hydroxylapatite Bone Grafting Material

Catalog Numbers	
Description	Catalog No.
Individual Pack, Includes 1.5Gr Ampule, 1 Small-Bore Syringe and 1 Funnel	0006
Individual Pack, Includes 4Gr Ampule, 1 Small-Bore Syringe and 1 Funnel	0007
Individual Pack, Includes 1.5Gr Pre-Filled Syringe	0008

Calcitite HA 4060 Particles

Catalog Numbers	
Description	Catalog No.
Four Sterile 0.5Gr Ampules	0003

Patient education made easy. And effective.

Zimmer Dental provides you with the tools essential to educate your patients about the benefits of dental implants. Patient education materials are a great way to begin the process.

Dental Implant Education

- Help patients accept dental implants as a treatment modality
- Allow the patient to be part of the decision process
- Explain why dental implants are not just an option, but represent the standard of care
- Help build trust in your treatment (can be built into your consultation)
- Help patients visualize the consequences of leaving a missing tooth untreated, such as bone resorption, tissue degeneration, or the structure of the mouth changing
- Provide office staff with the tools to participate in the presentation of dental implants

Contact your local Zimmer Dental Sales Consultant today to find out more about our Patient Education tools or consult our website at www.zimmerdental.com

Patient Education Brochures

Catalog Numbers		
Description	Catalog No.	Qty/Pkg
"Restoring Quality to Life" Brochure	0287	10/Pkg
"Ask Your Practitioner" Brochure Holder, Includes 1 Package of "Restoring Quality to Life" Brochures	4836	1 Ea
<i>Puros</i> Allografts "Healthy Bone Growth" Brochure	5590	25/Pkg

Patient Education Flip Chart

Catalog Numbers		
Description	Catalog No.	Qty/Pkg
Consider the Dental Implant Alternative	0293	1 Ea

Patient Education Models

Catalog Numbers		
Description	Catalog No.	Qty/Pkg
Treatment Planning Model (For Internal and External Dental Implant Connections)	2001	1 Ea
<i>SwissPlus</i> Implant and <i>Hader</i> Bar Overdenture Model	MSPHB	1 Ea
<i>Tapered Screw-Vent</i> ® Implant Maxillary Crown and Bridge Model	MTSVCB	1 Ea
<i>AdVent</i> ® Implant Mandibular Ball and Cap Overdenture Model	MAVOD	1 Ea
<i>Spline</i> ® Implant Deluxe Model	0871	1 Ea
<i>Locator</i> Implant Model (Internal Connection Implants)	LOCMDL	1 Ea

CAT NO.	DESCRIPTION	PAGE NO.	CAT NO.	DESCRIPTION	PAGE NO.
IMPLANTS					
TAPERED SWISSPLUS IMPLANTS WITH MTX SURFACE					
Measurements are listed as follows: Implant Body Diameter, Platform Diameter, Interface, Length.					
SPMB8	3.7mm, 3.8mm, 2.5mm Internal Hexagon, 8mm	2			
SPMB10	3.7mm, 3.8mm, 2.5mm Internal Hexagon, 10mm	2			
SPMB12	3.7mm, 3.8mm, 2.5mm Internal Hexagon, 12mm	2			
SPMB14	3.7mm, 3.8mm, 2.5mm Internal Hexagon, 14mm	2			
SPB8	3.7mm, 4.8mm, 3.0mm Internal Octagon, 8mm	2			
SPB10	3.7mm, 4.8mm, 3.0mm Internal Octagon, 10mm	2			
SPB12	3.7mm, 4.8mm, 3.0mm Internal Octagon, 12mm	2			
SPB14	3.7mm, 4.8mm, 3.0mm Internal Octagon, 14mm	2			
SPWB8	4.8mm, 4.8mm, 3.0mm Internal Octagon, 8mm	2			
SPWB10	4.8mm, 4.8mm, 3.0mm Internal Octagon, 10mm	2			
SPWB12	4.8mm, 4.8mm, 3.0mm Internal Octagon, 12mm	2			
SPWB14	4.8mm, 4.8mm, 3.0mm Internal Octagon, 14mm	2			
STRAIGHT SWISSPLUS IMPLANTS WITH MTX SURFACE					
Measurements are listed as follows: Implant Body Diameter, Platform Diameter, Interface, Length.					
OPB8	4.1mm, 4.8mm, 3.0mm Internal Octagon, 8mm	2			
OPB10	4.1mm, 4.8mm, 3.0mm Internal Octagon, 10mm	2			
OPB12	4.1mm, 4.8mm, 3.0mm Internal Octagon, 12mm	2			
OPB14	4.1mm, 4.8mm, 3.0mm Internal Octagon, 14mm	2			
OPWB8	4.8mm, 4.8mm, 3.0mm Internal Octagon, 8mm	2			
OPWB10	4.8mm, 4.8mm, 3.0mm Internal Octagon, 10mm	2			
OPWB12	4.8mm, 4.8mm, 3.0mm Internal Octagon, 12mm	2			
OPWB14	4.8mm, 4.8mm, 3.0mm Internal Octagon, 14mm	2			
SURGICAL COVER SCREWS					
Measurements are listed as follows: Platform Diameter, Component Length, Flare Diameter.					
SPSCS	3.8mm, 1.0mm, 4.6mm	2			
SPSC	3.8mm, 2.0mm, 4.6mm	2			
SPSCL	3.8mm, 3.0mm, 4.6mm	2			
OPSCS	4.8mm, 1.0mm, 5.2mm	2			
OPSC	4.8mm, 2.0mm, 5.2mm	2			
OPSCL	4.8mm, 3.0mm, 5.2mm	2			
OPWSC	4.8mm, 2.0mm, 6.0mm	2			
OPWSCS	4.8mm, 3.0mm, 6.0mm	2			
PROSTHETICS					
TRANSFER COMPONENTS					
OPT/5	Indirect Transfer, 4.8mm, 5.2mm	3			
OPT/6	Indirect Transfer, 4.8mm, 6.0mm	3			
SPMT	Indirect Transfer, 3.8mm, 4.6mm	3			
OPR	Replica w/Internal Octagon, 4.8mm	3			
SPMR	Replica w/Internal Hexagon, 3.8mm	3			
OPAR	Narrow Abutment/Shoulder Replica, 4.8mm	3			
SPMAR	Narrow Abutment/Shoulder Replica, 3.8mm	3			
STRAIGHT ABUTMENTS					
OPA	Narrow Straight Abutment w/Plastic Transfer/Waxing Coping, 4.8mm	4			
OPA2	Narrow Straight Abutment w/Plastic Transfer/Waxing Coping, 4.8mm	4			
SPMA2	Narrow Straight Abutment w/Plastic Transfer/Waxing Coping, 3.8mm	4			
OPC	Replacement Plastic Waxing Coping, 4.8mm	4			
SPMC	Replacement Plastic Waxing Coping, 3.8mm	4			
OPA/5	Straight Abutment/Transfer, Engages Internal Octagon, 4.8mm, 5.2mm	4			
OPA/6	Straight Abutment/Transfer, Engages Internal Octagon, 4.8mm, 6.0mm	4			
FMTM2	Straight Abutment/Transfer, Engages Internal Hexagon, 3.8mm, 4.6mm	4			
GPCAS	Replacement Screw for Straight Abutment, 4.8mm	4			
MHLAS	Replacement Screw for Narrow Straight Abutment, 3.8mm	4			
OPAS	Replacement Transfer Screw for Straight Abutment/Transfer, 4.8mm	4			
URS3	Replacement Transfer Screw for Straight Abutment/Transfer, 3.8mm	4			
FMS	Fixture Mount/Transfer Regular, 4.8mm, 5.2mm	5			
FMSW	Fixture Mount/Transfer Wide, 4.8mm, 6.0mm	5			
FMSM	Fixture Mount/Transfer Mini, 3.8mm, 4.6mm	5			
URS1	Fixture Mount Screw (Regular and Wide)	5			
URS2	Fixture Mount Screw (Mini)	5			
20° ANGLED ABUTMENTS					
OPH20	Internal Octagon, 4.8mm	5			
OPH20S	Replacement Screw for Internal Octagon	5			
SPH20	Internal Hexagon, 3.8mm	5			
AH20S	Replacement Screw for Internal Hexagon	5			
"CAST-TO" GOLD ABUTMENTS, ENGAGING					
OPGA	Gold Abutment, 4.8mm	6			
GPCAS	Replacement Screw, 4.8mm	6			
WSX	Waxing Screw, 4.8mm	6			
SPGA	Gold Abutment, 3.8mm	6			
MHLAS	Replacement Screw, 3.8mm	6			
MTWSD	Waxing Screw, 3.8mm	6			
OPS	Castable Plastic Sheath	6			
"CAST-TO" GOLD & PLASTIC COPINGS, NON-ENGAGING					
OPGC	Gold Coping, 4.8mm, 5.0mm	6			
OPCC	Plastic Castable Coping, 4.8mm, 5.0mm	6			
TAPERED ABUTMENTS					
OPACT	SwissPlus Tapered Abutment, 4.8mmD, 1.6mL	7			
OPACT3	SwissPlus Tapered Abutment, 4.8mmD, 3mL	7			
TAPERED HEALING CAP					
TATHC	Tapered Abutment Titanium Healing Cap	7			
TAPERED ABUTMENT TRANSFER COMPONENTS					
ACTDT	Direct Transfer	8			
ACTIT	Indirect Transfer	8			
ACTR	Abutment Replica	8			
TAPERED ABUTMENT COPINGS, SHEATHS & BAR SYSTEMS					
ACTGC	Abutment Collar Gold Coping	8			
ACTP	Plastic Castable Coping, 5.0mm	8			
TGC3	Bar Gold Coping, 3.0mm	8			
TGC5	Bar Gold Coping, 5.0mm	8			
TTC5	Titanium Bar Coping, 5.0mm	8			
ACTT	Titanium Temporary Coping, 9.5mm	8			
SCTS	Replacement Screw for Copings	8			
SCWS	Abutment Waxing Screw	8			
HGB	1.8mm Round Gold Bar, 50mm	8			
DGB	Dolder Resilient Gold Bar, 50mm	8			
BS1	Hader Bar Clip System	8			
LOCATOR ABUTMENTS FOR SWISSPLUS IMPLANTS					
Measurements are listed as follows: Implant Platform Diameter, Cuff Height.					
OPLOC4/1	4.8mm, 1.29mm	9			
OPLOC4/2	4.8mm, 2mm	9			
OPLOC4/3	4.8mm, 3mm	9			
OPLOC4/4	4.8mm, 4mm	9			
OPLOC3/1	3.8mm, 0.75mm	9			
OPLOC3/2	3.8mm, 2mm	9			
OPLOC3/3	3.8mm, 3mm	9			
OPLOC3/4	3.8mm, 4mm	9			
LOCATOR ATTACHMENT COMPONENTS					
LOCMP2	Male Processing Package, 2 Ea	9			
LOCMP10	Male Processing Package, 10 Ea	9			
LOCPP	Paralleling Post	9			
LOCRMW	Replacement Male	9			
LOCRMP	Replacement Male, Light Retention	9			
LOCRMB	Replacement Male, Extra Light Retention	9			
LOCRMG	Extended Range Replacement Male	9			
LOCRMRR	Extended Range Replacement Male, Light Retention	9			
LOCATOR ABUTMENT TRANSFER COMPONENTS					
LOCIMP	Impression Coping	10			
LOCAN/4	Female Analog, 4mmD	10			
LOCAN/5	Female Analog, 5mmD	10			
LOCATOR BAR DRILL AND TAP COMPONENTS					
LOCBF2	Locator Bar Female (2mm Thread)	10			
LOCBPP	Locator Yellow Bar Processing Package	10			
LOCDBT	Locator Bar Drill and Tap Holder	10			
LOCCTI	Locator 2.0 mm Castable Threaded Insert (Optional)	10			
LOCBT2	Locator 2.0mm Bar Tap (2mm Thread)	10			
LOCBD	Locator 1.7mm Bar Drill (2mm Thread)	10			

CAT NO.	DESCRIPTION	PAGE NO.	CAT NO.	DESCRIPTION	PAGE NO.
LOCATOR BAR LASER WELD COMPONENTS					
LOCLBF	Locator Laser Bar Female	10	10100802	S-8 O Oscillating Saw	13
LOCBPP	Locator Yellow Bar Processing Package	10	10100803	S-8 S Sagittal Saw	13
LOCPMAND	Locator Paralleling Mandrel	10	04726000	Surgical Saw Blades, Reciprocating, 15mm (1/Pkg)	13
			04800100	Surgical Saw Blades, Reciprocating, 15mm (5/Pkg)	13
			04725900	Surgical Saw Blades, Reciprocating, 20mm (1/Pkg)	13
			04800000	Surgical Saw Blades, Reciprocating, 20mm (5/Pkg)	13
			04725700	Surgical Saw Blades, Reciprocating, 30mm (1/Pkg)	13
			04799900	Surgical Saw Blades, Reciprocating, 30mm (5/Pkg)	13
			04800500	Surgical Saw Blades, Oscillating, 8mm (1/Pkg)	13
			04800400	Surgical Saw Blades, Oscillating, 15mm (1/Pkg)	13
			04726200	Surgical Saw Blades, Sagittal, 6mm (1/Pkg)	13
			04800300	Surgical Saw Blades, Sagittal, 6mm (5/Pkg)	13
			04726100	Surgical Saw Blades, Sagittal, 10mm (1/Pkg)	13
			04800200	Surgical Saw Blades, Sagittal, 10mm (5/Pkg)	13
LOCATOR BAR CAST-TO COMPONENTS					
LOCCTB	Locator Cast-To Bar Female	10			
LOCBPP	Locator Yellow Bar Processing Package	10			
LOCPMAND	Locator Paralleling Mandrel	10			
LOCATOR ABUTMENT TOOLS					
LOCCT2	Locator Core Tool	10			
LOCCTR	Locator Core Tool Replacement Tip (Fits LOCCT)	10			
LOCTW15	Locator Torque Wrench Insert Driver, 15mmL	10			
LOCTW21	Locator Torque Wrench Insert Driver, 21mmL	10			
BALL ABUTMENTS					
OPBA	Ball Abutment & Cap Attachment, 4.8mmD, 1.6mmL	11			
OPBA3	Ball Abutment & Cap Attachment, 4.8mmD, 3.0mmL	11			
BALL ABUTMENT TRANSFERS & REPLICAS					
OPBAT	Ball Abutment Transfer	11			
OPBAR	Ball Abutment Replica	11			
OVERDENTURE ATTACHMENT PROSTHETIC COMPONENTS					
CAS	Cap Attachment System	11			
CAI	Cap Attachment Instruments	11			
CA	Cap Attachment Housing & Cap Attachment	11			
CAH	Cap Attachment Housing	11			
CAN	Cap Attachment Nylon Liner (Transparent)	11			
CAN-G	Cap Attachment Nylon Liner (Grey-Rigid Retention)	11			
CAT	Cap Attachment Transfer	11			
CAB	Castable Ball Pattern	11			
CANM	Micro Cap Attachment Nylon Liner	11			
PROSTHETIC TOOLS AND TORQUE WRENCHES					
HX1.25	1.25mmD Hex Tool, 17mmL	12			
HXL1.25	1.25mmD Hex Tool, 22mmL	12			
THX1.25	1.25mmD Stainless Steel Tapered Frictional, 27mmL	12			
THXL1.25	1.25mmD Long Stainless Steel Tapered Frictional, 27mmL	12			
HX1.25D	Hex Drill	12			
TW30	Torque Wrench, 30 Ncm	12			
TW20	Torque Wrench, 20 Ncm	12			
TW1.25	Torque Wrench Hex Driver, for TW20 & TW30	12			
TW1.25L	Long Torque Wrench Hex Driver, for TW20 & TW30	12			
LOCTW15	Locator Torque Wrench Insert Driver, 15mmL, for TW20 and TW30	12			
LOCTW21	Locator Torque Wrench Driver, 21mmL, for TW20 and TW30	12			
LOCCT2	Locator Core Tool	12			
LOCCTR	Locator Core Tool Replacement Tip (Fits LOCCT)	12			
2320	Prosthetic Kit Tray, Autoclavable	12			
PTS20	PTS Disposable Inserts, 20 Ncm, White	12			
PTS30	PTS Disposable Inserts, 30 Ncm, Black	12			
ABTH	Abutment Holder	12			
PR	Reamer for Tapered Abutment Copings & SwissPlus Copings & Abutments	12			
SRT	Removal Tool for Abutment Screws or Fixation Screws	12			
	w/1.25mm Hex, Carbide Steel	12			
TLRT2	Removal Tool for Internal Hex Implant Abutments	12			
HLRTX2	Removal Tool for Internal Octagon	12			
SURGICAL ACCESSORIES					
ZIMMER SURGICAL MOTOR SYSTEM					
00900125	Zimmer Surgical Motor System (100-130V)	13			
00900124	Zimmer Surgical Motor System (220-240V)	13			
SURGICAL HANDPIECES AND SURGICAL SAW BLADES					
10207557	WI-75 E/KM Contra-Angle Handpiece 20:1	13			
10209201	WS-92 Contra-Angle Handpiece 1:2	13			
00001100	S-11 Straight Handpiece 1:1	13			
10101200	S-12 Angled Handpiece 1:2	13			
10100801	S-8 R Reciprocating Saw	13			
IRRIGATION ACCESSORIES					
04363600	One-Piece Irrigation Tubing Set, Disposable	14			
04719400	One-Piece Irrigation Tubing Set, Autoclave	14			
0285540R	Irrigation Spike with Roller, Autoclave	14			
04013900	Pump Tubing with Connectors, Autoclave	14			
04014000	Spare Pump Tubing without Connectors, Autoclave	14			
00929300	Irrigation Tubing, 2.2m, Disposable	14			
02610500	Internal Bur Coolant Tube	14			
02653901	Y-Connector	14			
02654700	Y-Tube Assembly	14			
04019000	Set of Tube Clamps	14			
04757100	Spray Clip for External/Internal Irrigation	14			
04005900	Irrigation Stand	14			
MICROMOTORS AND ACCESSORIES					
04720022	Motor with Cable for Surgical Motor System, 1.8m	14			
04009600	Motor with Cable for ImplantMED Surgical Motors, 1.8m	14			
04735700	Motor Holder for Zimmer Surgical Motor System	14			
02038200	Oil Spray Cap for Handpieces	14			
02139800	Testing Gauge for Burs	14			
04032600	Sterilizing Protector for Motor	14			
04013500	Sterilizing Box for Motor	14			
04013600	Transportation Case	14			
04541900	ST-GC Treatment Trolley with Shelf	14			
ASEPTIC SURGICAL DRAPE KIT					
ASD	Complete Surgical Drape Kit	14			
SURGICAL TOOLS					
SWISSPLUS IMPLANT SURGICAL TRAYS					
OPCST	Complete Surgical Kit	15			
OPSUR	Surgical Tray	15			
DRIVA DRILLS					
SV2.3DS	2.3mmD, 11mmL	18			
OP2.3D	2.3mmD, 17mmL	18			
SV2.8DS	2.8mmD, 11mmL	18			
OP2.8D	2.8mmD, 17mmL	18			
TLSDS	3.5mmD, 11mmL	18			
OP3.5D	3.5mmD, 17mmL	18			
SVWDS	4.2mmD, 11mmL	18			
OP4.2D	4.2mmD, 17mmL	18			
0201	Tapered Pilot Drill, 2.1mmD/1.6mmD, 8mmL	18			
DRIVA STEP DRILLS					
TSV3DS	3.4mm/2.8mmD, 11mmL	18			
SP3D	3.4mm/2.8mmD, 17mmL	18			
TSV4DS	4.4mm/3.8mmD, 11mmL	18			
SP4D	4.4mm/3.8mmD, 17mmL	18			
CLEANING TOOL					
NM1940	Drill Cleaning Wire	18			
CORTICAL BONE TAPS					
DT3.7	3.7mm Double-Lead Cortical Tap Tool	18			
DT4.8	4.8mm Double-Lead Cortical Tap Tool	18			

CAT NO.	DESCRIPTION	PAGE NO.	CAT NO.	DESCRIPTION	PAGE NO.
SINGLE PATIENT DRILLS			CALCITITE HA 2040 PARTICLES – NONRESORBABLE DENSE HYDROXYLAPATITE BONE GRAFTING MATERIAL		
DSSK	Drill Kit, Single Patient 3.7mm <i>Tapered SwissPlus</i> Implant.....	19	0006	Individual Pack, includes 1.5Gr Ampule, 1 Small-Bore Syringe and 1 Funnel.....	23
DSTK3S	Drill Kit, Single Patient 3.7mm <i>Tapered SwissPlus</i> Implant.....	19	0007	Individual Pack, includes 4Gr Ampule, 1 Small-Bore Syringe and 1 Funnel.....	23
DSSW	Drill Kit, Single Patient 4.8mm <i>Tapered SwissPlus</i> Implant.....	19	0008	Individual Pack, includes 1.5Gr Pre-Filled Syringe.....	23
DST4S	Drill Kit, Single Patient 4.8mm <i>Tapered SwissPlus</i> Implant.....	19	CALCITITE HA 2040 PARTICLES		
HEX TOOLS			0003	Four Sterile 0.5Gr Ampules.....	23
HX1.25	1.25mmD Hex Tool, 17mmL.....	20	EDUCATION PRODUCTS		
HXL1.25	1.25mmD Hex Tool, 22mmL.....	20	PATIENT EDUCATION BROCHURES		
THX1.25	1.25mmD Tapered Frictional Hex Tool, 23mmL.....	20	0287	“Restoring Quality to Life” Brochure.....	25
THXL1.25	1.25mmD Tapered Frictional Hex Tool, 27mmL.....	20	4836	“Ask Your Practitioner” Brochure Holder.....	25
HX1.25D	1.25mmD Hex Drill.....	20	5590	<i>Puros</i> Allografts “Healthy Bone Growth” Brochure.....	25
SURGICAL ARMAMENTARIA			PATIENT EDUCATION FLIP CHART		
SSHS	Stainless Steel Screwdriver Handle.....	20	0293	Consider the Dental Implant Alternative.....	25
RSR	<i>GemLock</i> Retaining Square Ratchet.....	20	PATIENT EDUCATION MODELS		
DE	Drill Extender.....	20	2001	Treatment Planning Model.....	25
PPAR	Paralleling Tool.....	20	MSPHB	<i>SwissPlus</i> Implant and <i>Hader</i> Bar Overdenture Model.....	25
OT3.0-S	3.0mmD Octagon Insertion Tool, Short, 17mmL.....	20	MTSVCB	<i>Tapered Screw-Vent</i> Implant Maxillary Crown and Bridge Model.....	25
OTL3.0-S	3.0mmD Octagon Insertion Tool, Long, 29mmL.....	20	MAVOD	<i>AdVent</i> Implant Mandibular Ball and Cap Overdenture Model.....	25
RHD2.5	2.5mmD <i>GemLock</i> Retaining Hex Drill, 21mmL.....	20	0871	<i>Spline</i> Implant Deluxe Model.....	25
RH2.5	2.5mmD <i>GemLock</i> Retaining Hex Driver, Short, 17mmL.....	20	LOCMDL	<i>Locator</i> Implant Model.....	25
RHL2.5	2.5mmD <i>GemLock</i> Retaining Hex Driver, Long, 28mmL.....	20	ZIMMER REGENERATIVE PRODUCTS		
IRT	Removal Tool for Screw-Type Implants, Carbide Steel.....	20	PUROS BONE ALLOGRAFTS		
HLRTX2	Removal Tool for Internal Octagon Implant Interface.....	20	8210R	<i>Puros</i> Cancellous Particles, 0.5cc, 250-1000.....	21
ZIMMER REGENERATIVE PRODUCTS			8211R	<i>Puros</i> Cancellous Particles, 1.0cc, 250-1000.....	21
PUROS BONE ALLOGRAFTS			8209R	<i>Puros</i> Cancellous Particles, 2.0cc, 250-1000.....	21
8210R	<i>Puros</i> Cancellous Particles, 0.5cc, 250-1000.....	21	8212R	<i>Puros</i> Cancellous Particles, 0.5cc, 1000-2000.....	21
8211R	<i>Puros</i> Cancellous Particles, 1.0cc, 250-1000.....	21	8213R	<i>Puros</i> Cancellous Particles, 1.0cc, 1000-2000.....	21
8209R	<i>Puros</i> Cancellous Particles, 2.0cc, 250-1000.....	21	8214R	<i>Puros</i> Cancellous Particles, 2.0cc, 1000-2000.....	21
8212R	<i>Puros</i> Cancellous Particles, 0.5cc, 1000-2000.....	21	8271R	<i>Puros</i> Cortical Particles, 0.5cc, 250-1000.....	21
8213R	<i>Puros</i> Cancellous Particles, 1.0cc, 1000-2000.....	21	8272R	<i>Puros</i> Cortical Particles, 1.0cc, 250-1000.....	21
8214R	<i>Puros</i> Cancellous Particles, 2.0cc, 1000-2000.....	21	8273R	<i>Puros</i> Cortical Particles, 2.0cc, 250-1000.....	21
8271R	<i>Puros</i> Cortical Particles, 0.5cc, 250-1000.....	21	8274R	<i>Puros</i> Cortical Particles, 0.5cc, 1000-2000.....	21
8272R	<i>Puros</i> Cortical Particles, 1.0cc, 250-1000.....	21	8275R	<i>Puros</i> Cortical Particles, 1.0cc, 1000-2000.....	21
8273R	<i>Puros</i> Cortical Particles, 2.0cc, 250-1000.....	21	8276R	<i>Puros</i> Cortical Particles, 2.0cc, 1000-2000.....	21
8274R	<i>Puros</i> Cortical Particles, 0.5cc, 1000-2000.....	21	8220R	<i>Puros</i> Block Allograft, 10mm.....	21
8275R	<i>Puros</i> Cortical Particles, 1.0cc, 1000-2000.....	21	8221R	<i>Puros</i> Block Allograft, 15mm.....	21
8276R	<i>Puros</i> Cortical Particles, 2.0cc, 1000-2000.....	21	PUROS MEMBRANE ALLOGRAFTS		
8220R	<i>Puros</i> Block Allograft, 10mm.....	21	8770R	<i>Puros</i> Pericardium Membrane, 15x20mm.....	22
8221R	<i>Puros</i> Block Allograft, 15mm.....	21	8771R	<i>Puros</i> Pericardium Membrane, 20x30mm.....	22
PUROS MEMBRANE ALLOGRAFTS			8772R	<i>Puros</i> Pericardium Membrane, 30x40mm.....	22
8770R	<i>Puros</i> Pericardium Membrane, 15x20mm.....	22	8793R	<i>Puros</i> Dermis Tissue Matrix, 10x10mm.....	22
8771R	<i>Puros</i> Pericardium Membrane, 20x30mm.....	22	8790R	<i>Puros</i> Dermis Tissue Matrix, 10x20mm.....	22
8772R	<i>Puros</i> Pericardium Membrane, 30x40mm.....	22	8791R	<i>Puros</i> Dermis Tissue Matrix, 10x40mm.....	22
8793R	<i>Puros</i> Dermis Tissue Matrix, 10x10mm.....	22	8792R	<i>Puros</i> Dermis Tissue Matrix, 20x40mm.....	22
8790R	<i>Puros</i> Dermis Tissue Matrix, 10x20mm.....	22	ZIMMER SOCKET REPAIR MEMBRANE		
8791R	<i>Puros</i> Dermis Tissue Matrix, 10x40mm.....	22	0154	<i>Zimmer</i> Socket Repair Membrane, 10x20 mm.....	22
8792R	<i>Puros</i> Dermis Tissue Matrix, 20x40mm.....	22	BIOMEND AND BIOMEND EXTEND ABSORBABLE COLLAGEN MEMBRANES		
ZIMMER SOCKET REPAIR MEMBRANE			0103	<i>BioMend</i> Membrane, 15x20mm.....	22
0154	<i>Zimmer</i> Socket Repair Membrane, 10x20 mm.....	22	0105	<i>BioMend</i> Membrane, 20x30mm.....	22
BIOMEND AND BIOMEND EXTEND ABSORBABLE COLLAGEN MEMBRANES			0107	<i>BioMend</i> Membrane, 30x40mm.....	22
0103	<i>BioMend</i> Membrane, 15x20mm.....	22	0140	<i>BioMend Extend</i> Membrane, 15x20mm.....	22
0105	<i>BioMend</i> Membrane, 20x30mm.....	22	0141	<i>BioMend Extend</i> Membrane, 20x30mm.....	22
0107	<i>BioMend</i> Membrane, 30x40mm.....	22	0142	<i>BioMend Extend</i> Membrane, 30x40mm.....	22
0140	<i>BioMend Extend</i> Membrane, 15x20mm.....	22	COLLATAPE, COLLACOTE AND COLLAPLUG COLLAGEN WOUND DRESSINGS		
0141	<i>BioMend Extend</i> Membrane, 20x30mm.....	22	0100	<i>CollaTape</i> Wound Dressing, 10/Pkg.....	23
0142	<i>BioMend Extend</i> Membrane, 30x40mm.....	22	0101	<i>CollaCote</i> Wound Dressing, 10/Pkg.....	23
COLLATAPE, COLLACOTE AND COLLAPLUG COLLAGEN WOUND DRESSINGS			0102	<i>CollaPlug</i> Wound Dressing, 10/Pkg.....	23
0100	<i>CollaTape</i> Wound Dressing, 10/Pkg.....	23			
0101	<i>CollaCote</i> Wound Dressing, 10/Pkg.....	23			
0102	<i>CollaPlug</i> Wound Dressing, 10/Pkg.....	23			

**For more information about our Products, Professional Programs
and Educational Opportunities, contact us:**

1900 Aston Avenue
Carlsbad, CA 92008-7308, USA
In the U.S. 800 854 7019

To fax an order 888 225 2483
Outside the U.S. +1 760 929 4300
Australia +61 (0)2 9950 5444 or 1 800 241 916
Canada +1 905 567 2073 or 1 800 265 0968

France +33 (0)1 45 12 35 35
Germany +49 (0)761 4584 722/723
Israel +972 (0)3 6124242
Spain +34 93 846 05 43

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsletterSignUp.aspx



www.zimmerdental.com