

BellaTek® Work Order Form



Account Number: _____
 Lab/Customer: _____
 Patient Reference: _____
 Prescribing clinician Zip: _____

Bill To: _____
 Ship To: _____
 Phone: _____

Design Review Needed?
 Yes No

Please fill in email below:

Tissue Displacement Options

Anatomical Tissue Displacement

3

Blanching Ok Moderate Tissue Displacement

2

No Blanching Minimal Tissue Displacement(Default)

1

Healing Abutment - black circle Margin - outer gray line
 Option three may require tissue adjustment and relief for placement.

Final Abutment Clearance

2.0 mm

- Default for **Posterior** : 2.0 mm
- Default for **Anterior** : 1.0 mm

Margin Design

■ **Shoulder (Default)**

■ **Chamfer**

■ **Feather**

Start Here With Tooth Number

Tooth #	Implant		Abutment Type		Margin Placement (Please select Supra or Sub when choosing margins)				Margin Design			Tissue Displacement Please select from 3 design options.			Final Abutment Clearance		
	Implant Type	Platform Diameter (mm)	Titanium*	Titanium Nitride (B3i Connections)	Zirconia (TSY/Connections)	<input type="checkbox"/> Apply same setting for all abutments				Shoulder (default)	Chamfer	Feather	3	2	1	Min	Max
						B/F	D	M	L								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								
						<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub	<input type="checkbox"/> Supra <input type="checkbox"/> Sub								

Parallel Abutments
 Connect blue dots for parallel abutments.

Spilted Restorations
 Connect blue dots for spilted restorations.

* If abutment type not selected, default to titanium.

Special Instructions: _____

Certification (by submitting this Work Order Form you certify the following):

- The stated information is correct, the submitted materials are accurate and do not contain metal.
- All items that have contacted the oral environment have been decontaminated.
- I have reviewed the applicable instructions (INST1079, INST1080, INST1120, INST1147, INST1156, INST1161) for this product.
- On behalf of the dentist (if not the dentist) that the soft-tissue has matured and healed completely.

This form authorizes the following:

- Fabrication of patient specific abutments
- Placement of analogs
- Modification of working models not consistent with applicable guidelines

All references to Zimmer Biomet Dental contained herein refer to the dental subsidiaries of Zimmer Biomet Holdings, Inc. Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured, distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). For additional product information, please refer to the individual product labeling/IFU. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet. ZB0024 REV A 11/16. ©2016 Zimmer Biomet, All rights reserved.